



Microsoft Update and PHP Update Monthly Notification

Summary: This notification contains information about Microsoft security updates released this month and approved for implementation on the ISM nStream and VaultStream platform. The only updates listed are those considered to be of a security or critical nature. Only the updates listed below have been validated with the nStream G3, nCareRX/MX, EasyLink Router, Thin Client, Status Board Engine, MedPresence Stand-Alone, All-In-One Client, and VaultStream systems.

Also included is information on the latest PHP updates that are approved for the v8.x and up VaultStream systems.

Product: nStreams, nCares, and Thin Clients connected to a network and VaultStream, as well as the Status Board Engine.

Date: August 15, 2018

Critical Security Bulletins (Obtained from Microsoft):

NOTE: Image Stream Medical **does not approve** MS updates **KB975560, KB975562, KB2631813 or KB2780091** and strongly suggests they be excluded from automatic update scenarios. Failure to exclude these updates will result in equipment performance issues.

Approved Updates:

August 2018 Updates		
KB Number	Platform	Title
nStream v5.4 Updates		
KB4343900	nStream	2018-08 Security Monthly Quality Rollup for Windows 7 for x86-based Systems
KB4345590	nStream	2018-08 Security and Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7
KB890830	nStream	Windows Malicious Software Removal Tool - August 2018
VaultStream v5.3, 8.0.x, v9.0.x Updates		
KB4343205	VaultStream	Cumulative Security Update for Internet Explorer 11 for Windows Server 2008 R2 for x64-based Systems
KB4343900	VaultStream	2018-08 Security Monthly Quality Rollup for Windows Server 2008 R2 for x64-based Systems
KB4345590	VaultStream	2018-08 Security and Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7 and Server 2008 R2 for x64
KB890830	VaultStream	Windows Malicious Software Removal Tool x64 - August 2018
EasyLink Router v8.0 Updates		
KB4339284	ELR v8.0	Update for WES09 and POSReady 2009
KB4340937	ELR v8.0	2018-08 Security Update for WES09 and POSReady 2009 for x86-based Systems
KB4343674	ELR v8.0	2018-08 Security Update for WES09 and POSReady 2009 for x86-based Systems
KB4344104	ELR v8.0	2018-08 Security Update for WES09 and POSReady 2009 for x86-based Systems
nCareRX v8.0 and v8.1, nCareMX/nStreamGX v8.0, v9.0, and Status Board Engine, 8.0 and v8.1 Updates		
KB4346744	nCare-SBE	2018-07 Update for .NET Framework 3.5.1 for Windows 7 and Server 2008 R2 for x64
KB4346407	nCare-SBE	2018-07 Update for .NET Framework 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7 and Server 2008 R2 for x64

KB4340556	nCare-SBE	2018-07 Security and Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7 and Server 2008 R2 for x64
KB4343900	nCare-SBE	2018-08 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x64-based Systems
	nCare-SBE	pciclearstalecache_1d566f78af356cc09c83cf173dc734c66659e0d2.exe
KB4343899	nCare-SBE	2018-08 Security Only Quality Update for Windows Embedded Standard 7 for x64-based Systems
KB4345590	nCare-SBE	2018-08 Security and Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7 and Server 2008 R2 for x64
kb4344152	nCare-SBE	windows6.1-kb4344152-x64_a4b728eaf3b80404339f0318e257e435758cb5fb.msu
KB4345679	nCare-SBE	2018-08 Security Only Update for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7 and Server 2008 R2 for x64
kb4344177	nCare-SBE	windows6.1-kb4344177-x64_ac4e288d33fead366bcbad5135d1aa6bc4d787f2.msu
HP Thin Client v8.0 and EasyLink Router v8.1 Updates		
KB4343899	ThinClient v8.0-ELR v8.1	2018-08 Security Only Quality Update for Windows Embedded Standard 7 for x86-based Systems
KB4343900	ThinClient v8.0-ELR v8.1	2018-08 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x86-based Systems
	ThinClient v8.0-ELR v8.1	pciclearstalecache_5d22303eb5f9d32d5d450051dcd7b7c6a55fdc2e.exe
MedPresence Stand-Alone v9.0.x and All-In-One Client (works with v9.0.x and 10.x)		
KB4091664	MedPresence-All-In-One	2018-07 Update for Windows 10 Version 1607 for x64-based Systems
KB4343887	MedPresence-All-In-One	2018-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4343902	MedPresence-All-In-One	2018-08 Security Update for Adobe Flash Player for Windows 10 Version 1607 for x64-based Systems
KB890830	MedPresence-All-In-One	Windows Malicious Software Removal Tool x64 - August 2018
KB4132216	MedPresence-All-In-One	2018-05 Update for Windows 10 Version 1607 for x64-based Systems
KB4103720	MedPresence-All-In-One	2018-05 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4103729	MedPresence-All-In-One	2018-05 Security Update for Adobe Flash Player for Windows 10 Version 1607 for x64-based Systems
KB4033631	MedPresence-All-In-One	Update for Windows 10 Version 1607 for x64-based Systems
KB4049411	MedPresence-All-In-One	Update for Windows 10 Version 1607 for x64-based Systems
	MedPresence-All-In-One	Intel Corporation - Bluetooth - 4/7/2017 12:00:00 AM - 19.60.0.3
ES4K v10.0.x - BlizCon		
KB4343902	BlizCon	2018-08 Security Update for Adobe Flash Player for Windows 10 Version 1607 for x64-based Systems
KB4343887	BlizCon	2018-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
ES4K v10.0.x - nCareRX2		
KB4343887	nCareRX2	2018-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4343902	nCareRX2	2018-08 Security Update for Adobe Flash Player for Windows 10 Version 1607 for x64-based Systems
ES4K v10.0.x - LSSB2		
KB4343887	LSSB2	2018-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems

KB4343902	LSSB2	2018-08 Security Update for Adobe Flash Player for Windows 10 Version 1607 for x64-based Systems
PHP Updates		
5.6.37 and 7.1.20		

For further questions concerning this notification please call Technical Support at (866) 376-6800.