

Microsoft® Update and PHP Update Monthly Notification

Summary: This notification contains information about Microsoft security updates released this month and approved for implementation on the Image Stream Medical EasySuite®, LiveStream®, nCare®, and VaultStream® platforms. The only updates listed are those considered to be of a security or critical nature. Only the updates listed below have been validated with the EasySuite and EasySuite 4K, nCare MX/RX, EasyLink Router, LiveStream StatusBoard Engine, and VaultStream systems.

Also included is information on the latest PHP updates that are approved for the v8.x and up VaultStream systems.

Products: EasySuites and nCares, connected to a network and VaultStream, as well as the LiveStream StatusBoard Engine.

Date: May 13, 2020

Critical Security Bulletins (Obtained from Microsoft):

NOTE: Image Stream Medical **does not approve** MS updates **KB975560, KB975562, KB2631813 or KB2780091** and strongly suggests they be excluded from automatic update scenarios. Failure to exclude these updates will result in equipment performance issues.

Approved Updates:

May 2020 Updates		
KB Number	Platform	Title
VaultStream v10.x Updates – Windows® Server 2012		
KB4550961	VaultStream 2012	2020-04 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems
KB4550958	VaultStream 2012	2020-04 Preview of Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems
KB890830	VaultStream 2012	Windows Malicious Software Removal Tool x64 - v5.82
KB4556401	VaultStream 2012	2020-05 Security and Quality Rollup for .NET Framework 3.5, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 8.1 and Server 2012 R2 for x64
KB4556846	VaultStream 2012	2020-05 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems
VaultStream v10.x Updates – Windows® Server 2016		
KB4550947	VaultStream 2016	2020-04 Cumulative Update for Windows Server 2016 for x64-based Systems
KB890830	VaultStream 2016	Windows Malicious Software Removal Tool x64 - v5.82
KB4556813	VaultStream 2016	2020-05 Cumulative Update for Windows Server 2016 for x64-based Systems
KB4552926	VaultStream 2016	2020-05 Cumulative Update for .NET Framework 4.8 for Windows Server 2016 for x64
nCare MX v10.2.2		
KB4552926	nCare MX	2020-05 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
KB4556813	nCare MX	2020-05 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
EasySuite 4K v10.3 EasySuite 4K Controller (a.k.a. BlizCon)		
KB4552926	EasySuite 4K Controller	2020-05 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
KB4556813	EasySuite 4K Controller	2020-05 Cumulative Update for Windows 10 Version 1607 for x64-based Systems

EasySuite 4K v10.3 EasySuite 4K Video Manager (a.k.a. nCare RX2)		
KB4552926	EasySuite 4K Video Manager	2020-05 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
KB4556813	EasySuite 4K Video Manager	2020-05 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
EasySuite 4K v10.3 – LiveStream StatusBoard 2		
KB4552926	LiveStream StatusBoard 2	2020-05 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
KB4556813	LiveStream StatusBoard 2	2020-05 Cumulative Update for windows 10 Version 1607 for x64-based Systems
MedPresence Standalone v10.3 (Frequently purchased with a custom cart)		
KB4552926	MedPresence Standalone v.10.3	2020-05 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
KB4556813	MedPresence Standalone v.10.3	2020-05 Cumulative Update for windows 10 Version 1607 for x64-based Systems
PHP Updates		
5.6.40 and 7.1.33 (7.3.18 available upon request for v10.2.2 VaultStreams only)		

For further questions concerning this notice, please call Technical Support at (866) 376-6800



Image Stream Medical

One Monarch Drive Littleton, MA 01460 U.S.A.
www.olympusamerica.com/olympus-systems-integration

Image Stream Medical, Inc., is a wholly owned subsidiary of Olympus Corporation of the Americas